INTENDED USE

The Surgis[®] Biodesign[™] Peyronie's Repair Graft is intended for implantation to reinforce soft tissue where weakness exists in the urological anatomy, including but not limited to urethroplasty, repair of tunica albuginea defects, and reinforcement in the repair of Peyronie's disease. The device is intended for one-time use.

[Rx ONLY] This symbol means the following:

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

This product is intended for use by trained medical professionals.

CONTRAINDICATIONS

This device is derived from a porcine source and should not be used for patients with known sensitivity to porcine materials.

PRECAUTIONS

- Do not resterilize. Discard all open and unused portions.
- Device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- Discard device if mishandling has caused possible damage or contamination, or if the device is past its expiration date.
- Ensure that device is rehydrated prior to suturing.
- Device performance has not been evaluated with suture spacing greater than 2 mm.
- Ensure that all layers of Peyronie's Repair Graft are secured when suturing.

POTENTIAL COMPLICATIONS

The following complications are possible with the use of surgical graft materials. If any of these conditions occur, the device should be removed:

- Infection
- Acute or chronic inflammation (Initial application of surgical graft materials may be associated with transient, mild, localized inflammation)
- Allergic reaction

STORAGE

This device should be stored in a clean, dry location at room temperature.

STERILIZATION

This device has been sterilized with ethylene oxide.

SUGGESTED INSTRUCTIONS FOR USING PEYRONIE'S REPAIR GRAFT

These recommendations are designed to serve only as a general guideline. They are not intended to supersede the institutional protocols or professional clinical judgment concerning patient care.

NOTE: Handle sheets using aseptic technique, minimizing contact with latex gloves.

REQUIRED MATERIALS

- A sterile dish (kidney dish or other bowl)
- Sterile forceps
- Rehydration fluid: at least 50 mL of room temperature sterile saline or sterile lactated Ringer's solution for each Peyronie's Repair Graft sheet.

1. Using aseptic technique, remove the Peyronie's Repair Graft inner pouch from its outer bag, and place the inner pouch in the sterile field.
2. Open the inner pouch carefully, and aseptically remove the Peyronie's Repair Graft sheet with the sterile forceps.
3. Place the Peyronie's Repair Graft sheet into the sterile dish in the sterile field. (Multiple Peyronie's Repair Graft sheets may be rehydrated simultaneously in the same dish.)
4. Add enough rehydration fluid to the dish to immerse the Peyronie's Repair Graft sheet.
5. Allow Peyronie's Repair Graft sheet to rehydrate for at least ten minutes.
6. Prepare the graft site using standard surgical techniques.
7. Using aseptic technique, trim the Peyronie's Repair Graft sheet to fit the site. User feedback suggests that the graft should be sized approximately 30% larger than the site.

Note: An alternative method is to cut the Peyronie's Repair Graft sheet to size prior to rehydration. If this method is selected, be sure to rehydrate the Peyronie's Repair Graft sheet prior to suturing it into place. See step 5.

8. Using aseptic technique, transfer the Peyronie's Repair Graft sheet to the graft site and suture into place using resorbable suture, while avoiding excess tension.

Note: Surgical experience indicates that suturing Peyronie's Repair Graft sheets with close tissue approximation produces better outcomes.

9. Complete the standard surgical procedure.
10. Discard any unused portions of the Peyronie's Repair Graft sheet.

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© COPYRIGHT COOK BIOTECH SEPTEMBER/2007
West Lafayette, IN 47906 U.S.A.
AUSTRALIA PTY. LTD.